



Cannula Selection

AEROX CANNULA SELECTION DATA SHEET

Step 1: Select Your Cannula

Standard Cannula

inexpensive and light weight, does not interfere with headsets.

Available in adult and child sizes.

Adult size part number: 4110-722

Child size part number: 4110-850



Oxysaver® Cannula (mustache type)

Our most popular cannula offers maximum oxygen conservation.

Available in Adult size only

Part number: 4110-701





Oxysaver® Cannula (pendent type)

For users that prefer a conserving cannula with lighter weight of a traditional cannula.

Part number: 4110-807



Step 2: Select Your Flow Meter



*Adjustable metering valve at the flowmeter: 4110-631-1 thru -6

*Flowmeter without metering valve: 4110-617-1 thru -6

Originators of High Duration Oxygen Systems



Duration (approximate) for aerox® Systems Using Oxysaver™ Cannula

	10,000 Ft. (MSL)					15,000 Ft. (MSL)					18,000 Ft. (MSL)				
Cylinder Size	A 180L 6cf	C 240L 9cf	D 400L 13cf	E-M 700L 13cf	F 1000L 33cf	A 180L 6cf	C 240L 9cf	D 400L 13cf	E-M 700L 13cf	F 1000L 33cf	A 180L 6cf	C 240L 9cf	D 400L 13cf	E-M 700L 13cf	F 1000L 33cf
USERS	Hour of use														
1	12	16	26.7	46.7	66.7	6.7	8.9	14.8	25.9	37	4.6	6.2	10.3	17.9	25.6
2	6	8	13.3	23.3	33.3	3.3	4.4	7.4	13	18.5	2.3	3.1	5.1	9	8.5
3	4	5.3	8.9	15.6	22.2	2.2	3	4.9	8.6	12.3	1.5	2.1	3.4	6	5.1
4	3	4	6.7	11.7	16.7	1.7	2.2	3.7	6.5	9.3	1.2	1.5	2.6	4.5	6.4
5	2.4	3.2	5.3	9.3	13.3	1.3	1.8	3	5.2	7.4	0.9	1.2	2.1	3.6	5.1
6	2	2.7	4.4	7.8	11.1	1.1	1.5	2.5	4.3	6.2	0.8	1	1.7	3	4.3

