



CYLINDER ASSEMBLY, CARBON FOR COMPOSITE FIBER-WRAPPED CYLINDERS



Specifications

Part Number: 4110-202

Dimensions:

Height: See Table

Width: See Table

Weight: See Table

Capacity: See Table

Thread: 3/4-16 UNF

Description: Carbon Composite Fiber Wrapped Aluminum Lined Cylinder per DOT-SP 10945.

Note: Cylinders are cleaned for oxygen service.

Maintenance:

Cylinder: 5 year Hydro-Static Test.
15 year service life.

Dimensional Table

Part Number	A (Height)	B (Diameter)	C (Weight)	D (Capacity) (Oxygen)
4110-202-1		DATA PENDING		
4110-202-2	19.2"	6.79"	6.7 Lbs Max	40 scf
4110-202-3		DATA PENDING		
4110-202-4		DATA PENDING		